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PRINCIPIO

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Feldman, A. T. and D. Wolfe (2014). Tissue Processing and Hematoxylin and Eosin Staining', in Histopathology

PRINCIPLE

The Papanicolaou test, which owes its name to Dr Georgios Papanicolaou, a pioneer in cytology and cancer, is performed to diagnose cervical cancer or cancer in other organs. It makes it possible to detect changes in cells which may be cancer precursors. The first step in this technique is the staining of nuclei with haematoxylin. The second step is the staining of cytoplasm with an orange solution, which stains the mature cells and the keratinised cells with different intensities. The so-called polychrome solution, which is a mixture of eosin, light green SF and Bismarck brown, is used in the third step. The differentiation of simple squamous epithelium is shown with the polychrome solution. The Haematoxylin and Eosin stain is the most commonly used stain in histology.

Haematoxylin stains are basically formed from metallic complexes of oxidised haematoxylin with divalent and trivalent ions. The mechanism of staining is carried out between the covalent bonds of the metal-haematoxylin compound and the anionic radicals of the tissue. The anionic groups of phosphoric acid in DNA will intervene if there is nuclear material in the reaction. The selectivity of the stain increases with the acid medium of the dye, since it hinders the staining of amphoteric elements, such as the amino acids of proteins. The nuclei appear blue to dark violet.

Progressive and regressive staining can be performed with haematoxylin. The structures of the nuclei show greater differentiation in regressive staining and are more easily seen.

INTENDED USE

Family of in vitro diagnostic medical devices intended for nuclear staining in histocytology.

They are qualitative stains that provide information on the physiological or pathological state of the tissue sample when the staining obtained is interpreted by a trained professional. To obtain a diagnosis, the result is utilized alongside other information such as the patient's medical history, physical status as well as other clinical and laboratory data obtained.

For trained and specialized personnel

REAGENTS AND COMPOSITION

Harris's Haematoxylin

1 x 500 mL	Ref. 99 22 32
1 x 1000 mL	Ref. 99 59 60
1 x 5000 mL	Ref. 99 12 97
1 x 1000 mL (STAINER)	Ref. 99 71 60

Composition

Haematoxylin C.I. No. 75 290	0.5%
Aluminium Sulphate	6%
Ethylene glycol	15%

Carazzi's Haematoxylin

1 x 1000 mL	Ref. 99 16 43
1 x 5000 mL	Ref. 99 41 75
1 x 1000 mL (STAINER)	Ref. 99 41 46

Composition

Haematoxylin C.I. No. 75 290	0.1%
Aluminium/Potassium Sulphate	6%
Glycerin	20%

Mayer's Haematoxylin

1 x 1000 mL	Ref. 99 77 50
1 x 5000 mL	Ref. 99 06 38

Composition

Haematoxylin C.I. No. 75 290	0.5%
Aluminium Ammonium Sulphate	5%
Glycerin	30%

Weigert's Haematoxylin A

1 x 500 mL	Ref. 99 38 01
------------	---------------

Composition

Haematoxylin C.I. No. 75 290	1 %
Ethanol 99	0.1%

For use in Masson's Trichrome Stain (Ref. 995700).

Weigert's Haematoxylin B

1 x 500 mL	Ref. 99 38 13
------------	---------------

Composition

Iron trichloride	2 %
HCl	< 0.01 %

For use in Masson's Trichrome Stain (Ref. 995700).

STORAGE AND STABILITY

Reagents which are stored at 10-35°C and protected from light are stable until the expiry date stated on the label. Containers must always be kept tightly closed.

Over time and temperatures below 10°C could lead to the formation of light precipitate in some reagents. This does not affect their functionality. The stain should be filtered before use. It is recommended to bring the reagent to room temperature (15-30°C), agitate and filter before using.

in Use stability should be determined by each user discretion.

MATERIAL NOT SUPPLIED

General-purpose pathology laboratory material.

Haematoxylin stain. (See QCA referencies available on the website)

Different concentrations of ethanol. (See QCA referencies available on the website)

Xylene or Xylene Substitute. (See QCA referencies available on the website)

Mounting medium

Microscope

Papanicolaou stains can be used with a manual staining method or with automatic staining equipment (see QCA stainers available on the website). It is recommended to follow the instructions for use of the automatic equipment supplied by the manufacturer.

PRECAUTIONS

All reagents must be handled with care by trained technicians. Safety instructions are on the product label. The safety data sheet should be consulted before using the reagents. Waste disposal must be carried out according to the local regulations in force.

Any serious incident that has occurred in relation to the device shall be reported to QCA and the competent authority of the member state in which the user is established.

SAMPLE

Cytological and histological samples of different origins: gynaecological, sputum, liquid biopsies.

Handling of the samples must be carried out in accordance with the established protocols in each laboratory for the preparation of samples for staining using state-of-the-art technology.

Handle the samples with care due to their potentially infectious nature.

QUALITY CONTROL

It is recommended to follow the quality control practices defined by the international organisms and periodic staining with routine quality control slide containing tissue processed in a similar manner to the test specimens should be performed.

PREPARATION OF THE WORKING REAGENT:

The stain is ready for use, for manual staining and for use in the QCA STAINER.

PROCEDURE

PAPANICOLAOU STAIN (MANUAL)

1. The smears which have been spray-fixed should be immersed in water for 5-10 minutes so that the carboxymethyl cellulose (CMC) that they contain is completely dissolved
2. Immerse in Harris's Haematoxylin for 1 minute
3. Rinse with water for 15 seconds
4. Rinse with ethanol-HCl (0.25%) for 10 seconds
5. Rinse with water for 15 seconds
6. Rinse with ammonia water (0.05% in NH₃) for 5 seconds
7. Rinse with water for 15 seconds
8. Immerse in ethanol 96° for 10 seconds (twice, using a different container each time)
9. Immerse in OG-6 stain for 3 minutes
10. Rinse in ethanol 96° for 15 seconds (twice, using a different container each time)
11. Immerse in EA 50, EA 36 or EA 65 stain for 3 minutes
12. Rinse in ethanol 96° for 15 seconds (twice, using a different container each time)
13. Immerse in absolute ethanol for 15 seconds (twice, using a different container each time)
14. Immerse in ethanol/xylol (50%) for 15 seconds
15. Immerse in xylol for 15 seconds
16. Immerse in xylol for 2 minutes. Mount the preparation

HAEMATOXYLIN AND EOSIN STAIN

1. Dewax the preparation
2. Rinse with distilled water for 1 minute after the last immersion in ethanol
3. Immerse in Haematoxylin for 8 minutes
4. Rinse with water for 1 minute
5. Rinse with ethanol-HCl (0.5%) for 13 seconds
6. Rinse with water for 15 seconds
7. Rinse with ammonia water (0.05% in NH₃) for 30 seconds
8. Rinse with water for 1 minute
9. Immerse in Eosin 1% for 30 seconds
10. Immerse in ethanol 96° for 6 seconds (twice, using a different container each time)
11. Immerse in absolute ethanol for 6 seconds (three times, using a different container each time)
12. Immerse in xylol/eucalyptol for 6 seconds (twice, using a different container each time)
13. Mount the preparation

MAYER'S HAEMATOXYLIN AND EOSIN STAIN

1. Dewax the preparation
2. Rinse with distilled water for 5 minutes after the last immersion in ethanol
3. Immerse in Mayer's Haematoxylin for 10 minutes
4. Rinse with water for 5 minutes
5. Immerse in alcoholic Eosin 1% for 2 minutes
6. Immerse in ethanol 96° for 10-15 seconds (twice, using a different container each time)
7. Immerse in absolute ethanol for 10-15 seconds (twice, using a different container each time)
8. Immerse in xylol (three times, using a different container each time)
9. Mount the preparation

GB - HARRIS'S HEMATOXYLIN

—

GB - CARAZZI'S HEMATOXYLIN

—

GB - MAYER'S HAEMATOXYLIN

—

For the QCA STAINER, the different methods are found in the system software.

Each user must validate this procedure in his laboratory in order to adjust it to his own standard method, applying the various variants of this procedure, both manual and automatic, allowing him to optimize the results. It must be taken into account that the intensity of the staining is proportional to the staining time.

RESULTS

Nuclei: blue-violet

Acidophilic cytoplasm: reddish

Basophilic cytoplasm: blue

Keratinised cytoplasm: orange

For Masson's Trichrome Stain see Ref. 995700

NOTES

The result of the staining should be treated as a guide and must be confirmed with additional tests.

The colour results indicated should be treated as a guideline as the intensity and shade of the colour may be influenced by several factors, such as the fixation and the staining time (the staining intensity is proportional to the staining time).

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